FOSFATASA ALCALINA P-NPP

MÉTODO DGKC

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Szasz, G., Rautenburg, H.W. (1971). Z. Kinderheilk., 111, 233 - 239.

Hausamen, T.U., Helger, R., Gross, W. (1967). Clin. Chim. Acta, 15, 241 - 245.

Recommendations of the German Society for Clinical Chemistry. (1970). Z. Klin. Chem. u Klin. Biochem., 8, 658 -660.

Tietz, NW., Textbook of Clinical Chemistry 6th Edition, W.B. Saunders, Philadelphia (2018).

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).


ALKALINE PHOSPHATASE P-NPP

DGKC METHOD

For “in vitro” determination of alkaline phosphatase in serum or plasma

PRINCIPLE

Determination of p-Nitrophenol produced by the interaction of Alkaline Phosphatase with the substrate p – Nitrophenylphosphate allows the measurement of enzymatic activity. In optimum conditions the $\Delta\text{Abs}/\text{min}$ value is directly related to the Alkaline Phosphatase concentration of the sample.

DIAGNOSTIC USE

Alkaline phosphatase is an enzyme whose concentration in serum may be altered by various diseases.

Its main clinical application is related to cases of obstructive liver disease, with significant elevations of basal level, especially in extrahepatic obstruction.

In metabolic bone diseases the increase of alkaline phosphatase level is associated to increased osteoblastic activity.

In bone metastasis the increase appears only in the case of osteosclerotic alterations.

Single test result can not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 1 x 125 mL. (Ref. 99 05 16). Contents:

A. 1 x 100 mL Buffer solution

Ref. 99 73 12

B. 1 x 25 mL Substrate

Ref. 99 68 58

WORKING REAGENT PREPARATION

Reagents A and B are ready-to-use. If a Monoreagent procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of A (Buffer solution) + 1 part of B (Substrate).

WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:

Diethanolamine buffer pH 9.8 1.2 M

p-nitrophenylphosphate 12 mM

Preservatives and stabilizers

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.

The Monoreagent is stable for 4 weeks at 2-8°C and for 3 week at room temperature ($\leq 25^\circ\text{C}$), when protected from the sunlight.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank ≥ 1.500 .

ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer thermostatable at 37°C. Cuvette: 1 cm light-path.

SAMPLE

Serum or heparinized plasma. Use samples free from hemolysis. Sera kept at 2-8°C will remain stable for 7 days.

CAUTION

The reagents contain sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the SDS before using the reagent.

The calibrator must be considered as a human sample, and, thus, potentially infectious. Use adequate protection.

The disposal of the residues has to be made according to local regulations.

QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

Bring the reagents and the analyzer to 37°C.

Monoreagent technique	37°C
Working reagent	1.0 mL
Sample	0.02 mL
Bireagent technique	37°C
Buffer solution (A)	1.0 mL
Sample	0.02 mL
Mix, incubate for approx. 1 minute	
Substrate (B)	0.25 mL

Mix, transfer to the reading cuvette and start the stopwatch. Read the absorbance after 1, 2, 3 and 4 min.

Reading

Wavelength: 405 nm

Blank: Water

Cuvette: Thermostated 1 cm light-path

CALCULATIONS

The formula indicated is used to obtain the factor to calculate the U/L

$$\Delta\text{Abs}/\text{min} \times \frac{V_t \times 10^6}{C \times l \times V_s} = \text{U/L}$$

Where:

Vt: Total volume of reaction

Vs: Sample volume

l: Cuvette light path

C: Extinction Coefficient of p-nitrophenol in alkaline medium, at 405 nm: 18500

Monoreagent 405 nm Factor: 2760

Bireagent 405 nm Factor: 3432

U/L = $\Delta\text{Abs}/\text{min} \times \text{Factor}$

REFERENCE VALUES

	Adults	Children 3-15 years
37°C	98-279 U/L	250-775 U/L

The stated values are for guidance. Each particular laboratory should establish its own normal range, using its own instrumentation, blood collection methods and test procedures.

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol.

These results have been obtained using a manual method.

Sensitivity, as detection limit: 6 U/L

Linearity: Up to 1000 U/L. For higher values, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.

Accuracy: 102%

Repetitivity, as Coefficient of Variation: 0.98%

Reproducibility, as Coefficient of Variation: 3.26%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

INTERFERENCES

Bilirubin (up to 30 mg/dL), intralipid (up to 5000 mg/dL), hemoglobin (up to 1000mg/dL), do not interfere.

Anticoagulants such as EDTA, oxalate or citrate which chelate divalent cations should not be used since they would result in enzyme inhibition.

REFERENCES

Szasz, G., Rautenburg, H.W. (1971). Z. Kinderheilk., 111, 233 - 239.

Hausamen, T.U., Helger, R., Gross, W. (1967). Clin. Chim. Acta, 15, 241 - 245.

Recommendations of the German Society for Clinical Chemistry. (1970). Z. Klin. Chem. u. Klin. Biochem., 8, 658 -660.

Tietz, N.W., Textbook of Clinical Chemistry 6th Edition, W.B. Saunders, Philadelphia (2018).

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

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GB - FOSF. ALCALINE LIQ. (DGKC), BUFFER

Danger

Hazard: Causes skin irritation. Causes serious eye damage. May cause damage to organs <or state all organs affected, if known> through prolonged or repeated exposure <state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard>.

Precautionary: Do not breathe dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor/... Specific treatment (see ... on this label). Take off contaminated clothing and wash it before reuse. Dispose of contents/container to ...

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