GAMMA - GT LÍQUIDA

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
GAMMA - GT LIQUID

SZASZ MODIFIED METHOD

For "in vitro" determination of gamma-GT in serum or plasma

PRINCIPLE

Serum γ -glutamyl transferase (γ -GT) transfer the glutamyl group to cosubstrate glycyl-glycine producing 5-amino-2-nitrobenzoate. In optimum conditions enzymatic activity is directly related to amino-2-nitrobenzoate produced that is measured by spectrophotometry.

DIAGNOSTIC USE

In 90% of liver diseases there is a rise of γ -GT levels. The values are moderate in cirrhosis and infectious mononucleosis, while significant increases were found in obstructive liver disease, intrahepatic cholestasis or primary biliary cirrhosis.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 1 x 50 mL. Ref. 99 35 61. Contents:

- A. 1 x 40 mL Buffer solution Ref. 99 92 48
B. 1 x 10 mL Liquid substrate Ref. 99 38 43

Kit 1 x 125 mL. Ref. 99 10 56. Contents:

- A. 1 x 100 mL Buffer solution Ref. 99 10 68
B. 1 x 25 mL Liquid substrate Ref. 99 10 81

Kit 1 x 250 mL. Ref. 99 10 54. Contents:

- A. 2 x 100 mL Buffer solution. Ref. 99 10 68
B. 1 x 50 mL Liquid substrate Ref. 99 10 82

WORKING REAGENT PREPARATION

Reagents A and B are ready-to-use. If a monoreagent procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of A (Buffer solution) + 1 part of B (Liquid substrate).

WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:

Tris buffer pH 8.4	90 mM
L-gamma-glutamyl-3-carboxy-p-nitroanilide	6 mM
Glycyl-glycine	100 mM
Stabilizers and preservatives	

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.

The monoreagent is stable 5 weeks at 2-8°C and 3 weeks at room temperature ($\leq 25^\circ\text{C}$), when protected from the sunlight.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank ≥ 1.75 .

ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer thermostatable at 37°C. Cuvette: 1 cm light-path.

CAUTION

The reagents contain sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the SDS before using the reagent. The calibrator must be considered as a human sample, and, thus, potentially infectious. Use adequate protection.

The disposal of the residues has to be made according to local regulations.

SAMPLE

Serum or plasma. Samples free from hemolysis should be used. The gamma-GT-asic activity will remain stable for 10 days if the sample is stored at 2-8°C.

QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

Bring the working reagent and the analyzer to the working temperature (25°C, 30°C, 37°C).

Monoreagent technique	25°C/30°C/37°C
Working reagent	1.0 mL
Sample	0.1 mL
Bireagent technique	25°C/30°C/37°C
Buffer sol.(A)	1.0 mL
Sample	0.1 mL
Mix, incubate for approx. 1 minute	
Substrate (B)	0.25 mL

Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2 and 3 min.

Determine the mean of Abs/min of different lectures.

Reading

Wavelength: 405 nm

Blank: Water

Cuvette: Thermostated 1 cm light path

CALCULATIONS

The formula indicated is used to obtain the factor to calculate the U/L

$$\Delta\text{Abs}/\text{min} \times \frac{V_t \times 10^6}{\epsilon \times l \times V_s} = \text{U/L}$$

V_t : Total reaction volume

V_s : Sample volume

l : Cuvette light path

ϵ : Extinction coefficient of 5-amino-2-nitrobenzoate at 405 nm: 9500

For the calculation of the U/L determine the average Abs/min and multiply by the factor:

Monoreagent method	U/L = $(\Delta\text{Abs}_{405\text{ nm}} / \text{min.}) \times 1158$
Bireagent method	U/L = $(\Delta\text{Abs}_{405\text{ nm}} / \text{min.}) \times 1420$

REFERENCE VALUES

Assay temp.	Men	Women
25 °C	6-28 U/L	4-18 U/L
30 °C	8-38 U/L	5-25 U/L
37 °C	11-50 U/L	7-32 U/L

The stated values are for guidance. Each particular laboratory should establish its own normal range, using its own instrumentation, blood collection methods and test procedures.

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using a manual method.

Sensitivity, as detection limit: 4 U/mL

Linearity: Up to 600 U/L. For higher values, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.

Accuracy: 98.4%

Repetitivity, as Coefficient of Variation: 1.77%

Reproducibility, as Coefficient of Variation: 2.29%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

INTERFERENCES

Hemolyzed samples will give false results. Fluoride, EDTA, citrate and oxalate inhibit the enzyme activity.

REFERENCES

Szasz, G., (1969) Clin. Chem., 15, 124 - 136.

Tietz, NW., Textbook of Clinical Chemistry 6th Edition, W.B. Saunders, Philadelphia (2018).

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