GOT / AST U.V. LÍQUIDA

MÉTODO IFCC

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GOT / AST U.V. LIQUID


IFCC METHOD

For “in vitro” determination of Aspartate aminotransferase GOT/AST in serum or plasma

PRINCIPLE

The glutamic-oxalacetic transaminase (GOT), actually AST, catalyzes the reaction between L-aspartic acid and α-ketoglutaric acid. The oxalacetic acid formed is reduced by NADH with the aid of auxiliary enzyme MDH. The change of NADH to NAD⁺ produces a change of Abs. The reaction media contains, also, LDH to remove endogenous pyruvate to avoid possible interferences. In optimum reaction conditions the ΔAbs/min is directly related to GOT concentration in the sample.

DIAGNOSTIC USE

Increases in GOT activity were observed in cases of liver damage: hepatitis of various types, necrosis or damage in hepatocytes, cholestatic ictericia. High levels are also seen in diseases affecting the heart muscle. In alcoholic hepatitis and acute myocardial infarction, the increase of GOT activity is greater than the activity of GPT.

Single test result can not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 1 x 50 mL. (Ref. 99 80 03). Contents:
A. 1 x 40 mL Enzymes solution Ref. 99 61 07
B. 1 x 10 mL Liquid substrate Ref. 99 22 00

Kit 1 x 250 mL. (Ref. 99 95 00). Contents:
A. 2 x 100 mL Enzymes solution Ref. 99 95 20
B. 1 x 50 mL Liquid substrate Ref. 99 21 65

Kit 1 x 940 mL. (Ref. 99 04 06). Contents:
A. 3 x 250 mL Enzymes solution Ref. 99 04 02
B. 1 x 190 mL Liquid substrate Ref. 99 04 11

WORKING REAGENT PREPARATION

Reagents A and B are ready-to-use. If a monoreagent procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of A (enzymes solution) + 1 part of B (liquid substrate).

WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:
Tris-HCl buffer pH 7.8 80 mM
L-aspartic acid 240 mM
α-ketoglutaric acid 12 mM
NADH 0.18 mM
MDH ≥ 600 U/L
LDH ≥ 800 U/L
Stabilizers and preservatives

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The Monoreagent is stable for 4 weeks at 2-8°C and for 1 week at room temperature (≤ 25°C), when protected from the sunlight.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank ≤ 1.0

ADDITIONAL EQUIPMENT

General laboratory equipment.
Spectrophotometer or photometer thermostatable at 37°C. Cuvette: 1 cm light-path.

SAMPLE

Serum or plasma with EDTA or heparin. Samples free from hemolysis should be used. Sera kept at 2-8°C loses ca. 10% of its activity after 3 days.

CAUTION

The reagents contain sodium azide at 0.09%. Handle with care.
The safety statements are on the label. It is advisable to look at the SDS before using the reagent.
The calibrator must be considered as a human sample, and, thus, potentially infectious. Use adequate protection.
The disposal of the residues has to be made according to local regulations.

QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 994148/996571) and Seriscann Anormal (Ref. 994685/999329) should be included in each test series. It is recommended to calibrate with the Calibrator for Autoanalyzers (Ref. 99 62 80) whenever the batch of reagent and/or calibrator is changed

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

The described method is the one proposed by the International Federation of Clinical Chemistry (IFCC). Bring reagents and the analyzer to working temperature.

Monoreagent technique	25/30°C	37°C
Working reagent	1.0 mL	1.0 mL
Sample	0.2 mL	0.1 mL
Bireagent technique	25/30°C	37°C
Enzymes sol.(A)	1.0 mL	1.0 mL
Sample	0.2 mL	0.1 mL
Mix, incubate for approx. 1 min		
Substrate (B)	0.25 mL	0.25 mL

Mix and start the stopwatch. Read the absorbance after 1, 2 and 3 min
Determine the mean of Abs/min of different lectures.

Reading

Wavelength: 334 nm; 340 nm; 365 nm
Blank: Water.
Cuvette: Thermostated, 1 cm light-path

CALCULATIONS

The formula indicated is used to obtain the factor to calculate the U/L

$$\Delta\text{Abs/min} \times \frac{V_t \times 10^6}{\epsilon \times l \times V_s} = \text{U/L}$$

Vt: Total reaction volume
Vs: Sample volume
l: Cuvette light path
ε: Extinction coefficient of NADH:
365 nm: 3.40 x 10³
340 nm: 6.31 x 10³
334 nm: 6.17 x 10³

Monoreagent procedure

	25/30°C	37°C
334 nm	ΔAbs/min x 970=U/L	ΔAbs/min x 1780=U/L
340 nm	ΔAbs/min x 950=U/L	ΔAbs/min x 1745=U/L
365 nm	ΔAbs/min x 1765=U/L	ΔAbs/min x 3235=U/L

Bireagent procedure

	25/30°C	37°C
334 nm	ΔAbs/min x 1175=U/L	ΔAbs/min x 2185=U/L
340 nm	ΔAbs/min x 1150=U/L	ΔAbs/min x 2140=U/L
365 nm	ΔAbs/min x 2130=U/L	ΔAbs/min x 3970=U/L

REFERENCE VALUES

Temperature	Men	Women
25 °C	≤ 18 U/L	≤ 15 U/L
30 °C	≤ 25 U/L	≤ 21 U/L
37 °C	≤ 37 U/L	≤ 31 U/L

The stated values are for guidance. Each particular laboratory should establish its own normal range, using its own instrumentation, blood collection methods and test procedures

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. The indicated results have been obtained using a automatic method for the BT 3500 autoanalyzer.

Sensitivity, as detection limit: 1.5 U/ml
Linearity: Up to 520 U/L. For higher values, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.
Accuracy: 96 %
Repetitivity, as Coefficient of Variation: 1.42%
Reproducibility, as Coefficient of Variation: 2.83%
Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.
Details of the performance studies are available on request

INTERFERENCES

Highly haemolysed sera will interfere with the reaction (from 500mg/dL). No bilirubin interference is observed up to 30mg/dL. When assaying high activity samples, a very low initial absorbance can be found which is mainly due to the rapid conversion of the NADH at the early stage of the reaction. In such a case, dilute the sample with saline (NaCl 0.9%) 1/10 and assay it once again. Multiply the final result by 10.

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